

Rapid detection of tetracycline resistance in bovine *Pasteurella multocida* isolates by MALDI Bityper antibiotic susceptibility test rapid assay (MBT-ASTRA)

Laura Van Driessche

Department of  
National Reference Laboratory for  
Diseases of Cattle

OPEN

Received: 9 May 2018

Accepted: 21 August 2018

Published online: 11 September 2018

SCIENTIFIC REPORTS



**Rapid detection of tetracycline resistance in bovine *Pasteurella multocida* isolates by MALDI Bityper antibiotic susceptibility test rapid assay (MBT-ASTRA)**


Laura Van Driessche<sup>1</sup>, Jade Bokma<sup>1</sup>, Linde Gille<sup>1</sup>, Pieter-Jan Ceyssens<sup>2</sup>, Katrin Sparbier<sup>3</sup>, Freddy Haesebrouck<sup>1,4</sup>, Piet Deprez<sup>1</sup>, Filip Boyen<sup>4</sup> & Bart Pardon<sup>1</sup>

<https://www.nature.com/articles/s41598-018-31562-8>

[Laura.vandriessche@ugent.be](mailto:Laura.vandriessche@ugent.be)

# INTRODUCTION

- Infectious bronchopneumonia  
leading cause antimicrobial  
use in calves
- ↗ Pressure intensive  
antimicrobial consumption food  
animals




# INTRODUCTION

- Reduction and rational antimicrobial use  
=> TOP PRIORITY!
- Formularies
  - 1<sup>st</sup> choice products
  - Certain classes: identification and susceptibility test requested




# INTRODUCTION






MALDI-TOF MS: Matrix-assisted laser desorption/ionization–time-of-flight mass spectrometry



or




‘fingerprint’

# OBJECTIVE

- Design and validate MBT-ASTRA for *Pasteurella multocida* (ubiquitous respiratory pathogen in cattle) and tetracycline (frequently used antimicrobial)
- Determine diagnostic accuracy of MBT-ASTRA and disk diffusion with MIC-gradient strip test



- MBT-ASTRA: Maldi Biotype-antibiotic susceptibility test rapid assay  
=> semi-quantitative MALDI-TOF MS to analyse susceptibility/ resistance status of bacteria in few hours





# DESIGN MBT-ASTRA




Conditions to determine:

1. Growth medium: **Cation-adjusted Mueller Hinton Broth**
2. Starting concentration of *P. multocida*:  **$10^7$  cfu/mL**
3. Concentration of tetracycline:  **$4\mu\text{g}/\text{mL}$**
4. Incubation time: **3 hours!!**






## Area under the curve (AUC)




## Relative growth ratio (RG)

$$\frac{\text{RG: AUC} + \text{Ab}}{\text{AUC} - \text{Ab}}$$




# DIAGNOSTIC ACCURACY MBT-ASTRA

- Database 100 recent clinical strains of *P. multocida* with MIC- values for tetracycline by MIC-gradient strip test



# DIAGNOSTIC ACCURACY MBT-ASTRA



# DIAGNOSTIC ACCURACY MBT-ASTRA AND DISK DIFFUSION

Test	Reference test (MIC-gradient strip test)	
	<i>Resistant</i>	<i>Susceptible</i>
MBT-ASTRA	<i>Resistant</i>	95.7% (44/46)
	<i>Susceptible</i>	4.3% (2/46)
	<i>Total</i>	100% (46/46)
Disk diffusion	<i>Resistant</i>	93.5% (43/46)
	<i>Susceptible</i>	6.5% (3/46)
	<i>Total</i>	100% (46/46)




# DIAGNOSTIC ACCURACY MBT-ASTRA AND DISK DIFFUSION



	Disk diffusion	MBT-ASTRA
<i>Essential agreement</i>	95%	98%
<i>Very major error</i>	3%	2%
<i>Major error</i>	2%	0%
<i>Minor error</i>	0%	0%
<i>Sensitivity</i>	93.5% (86.3%, 100.6%)	95.7% (89.8%, 101.5%)
<i>Specificity</i>	96.3% (91.3%, 101.3%)	100% (100.0%, 100.0%)
<i>RPV</i>	95.6% (89.5%, 101.6%)	100% (100.0%, 100.0%)
<i>SPV</i>	94.5% (88.5%, 100.5%)	96.4% (91.6%, 101.3%)

# CONCLUSION




**Identification and susceptibility testing possible  
in the same day!!**

# ACKNOWLEDGEMENTS



# CONCLUSION



<https://www.nature.com/articles/s41598-018-31562-8>